
REVOGYM

RP20

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(EN) SAFETY INSTRUCTIONS

1. Children may use the device only under adult supervision. Pregnant women and elderly persons should use the device with caution and in accordance with a doctor's recommendations.
2. Select the appropriate massage head depending on the body part being massaged.
3. Press and move the device gently. The massage time for the same area of the body should not exceed 1 minute.
4. High vibration speed should be used only on soft body tissues. Use caution. Do not use high speed or apply strong pressure on the head or bones.
5. If pain or discomfort is felt during the massage, stop using the device immediately.
6. Avoid direct contact with the protruding shaft and the heat-dissipation opening of the massage gun to prevent possible damage.
7. Do not disassemble the device
8. Use only the original charger and the included accessories.
9. If the device makes unusual sounds or the massage head becomes loose, switch off the device immediately.
10. To clean the surface of the device, use a slightly damp cloth, then dry it with a soft, dry cloth.
11. After each massage, remove the massage head and place it back in the case.

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(EN) USER MANUAL

Assembly Instructions

1. Select the appropriate massage head and insert it into the round opening at the front of the massager.
2. Turn on the main power switch located at the bottom of the handle.
3. Press the On/Off button once, located on the LCD screen.
4. Continue pressing the switch button to set the desired speed level.
5. After finishing the massage, remember to turn off the main power switch at the bottom of the handle.

Charging the Device

- Before first use, fully charge the device (approx. 3–4 hours).
- Do not use the device while charging to avoid potential battery damage.
- It is recommended to charge the device before the battery level drops below 20% to extend battery life.
- If the device is not used for a long time, the built-in battery should be charged regularly (about once a month).
- Do not charge the device in environments with high temperature or humidity, and do not expose it to direct sunlight.

Important Features

- **Timer Function:** If the massage gun is used continuously for 10 minutes, the device will automatically turn off to prevent motor overheating and battery damage.
- **12 Massage Heads & 7 Speed Levels:** The massage gun offers 12 interchangeable heads and 7 speed levels, allowing precise adjustment of massage intensity. Suitable for both low-intensity warm-up and deep muscle massage.

Frequently Asked Questions

1. Why won't the massage gun turn on?
 - Make sure the On/Off switch at the bottom of the handle is in the ON position.
 - The device may not start due to a low battery level. Try recharging it.
 - If the device still does not turn on, contact the service center or seller.
2. Why does the device become louder after some time of use?
 - During long or frequent use, the sound may increase, which is normal.
 - If there are noticeable unusual noises, loose screws, or damaged components, contact the service center or seller.
3. What is the function of the lamp on the rear cover?
 - The lamp indicates pressure: it does not light up when pressure is insufficient, and additional lights turn on as pressure increases.
 - The light may also indicate excessive pressure.
4. Is the device intelligent and can it increase power?
 - This function automatically increases rotational speed when resistance rises, ensuring more stable and powerful operation of the massager.

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