
REVOGYM

TP02

velaco.pl

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(EN) SAFETY INSTRUCTIONS

1. Use the device according to the manual. Do not allow children or unauthorized persons to operate it.
2. Operate only on a clean, dry, stable surface, away from moisture, water, flammable liquids, and fire sources.
3. Keep hair, clothing, and jewelry away from moving parts.
4. Always make sure the device is turned off before cleaning, changing attachments, or moving it.
5. Do not expose the device to water and do not block ventilation openings.
6. Regularly check the condition of the device and battery. Do not use if damaged.
7. Use only the supplied charger.
8. Maintain proper posture and balance during use.
9. Do not press too hard – use gently, max. 1 minute per body area.
10. After massage, wear dry clothing or dry skin with a towel.
11. Never leave the device on unattended.
12. For home use only.

Contraindications

Do not use if:

- you are pregnant,
- you suffer from diabetes, neuropathy, epilepsy, migraines, retinal problems, joint or spine disorders, or have metal implants,
- you have recently had surgery or joint replacement,
- you are under the influence of alcohol, drugs, or medication affecting control.

If in doubt, consult your physician.

Disposal Information

Our products are recyclable. At the end of their service life, please take the device to an appropriate waste disposal facility (local collection point).

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- (EN)
1. speed indicator
 2. battery indicator
 3. increase speed
 4. decrease speed
 5. touch switch
 6. power on/off switch
 7. charging ports
- Powerful brushless motor with high torque and quiet operation technology.
 - Replaceable battery suitable for recycling.
 - Lithium-ion battery.
 - 30-speed settings for optimal pressure adjustment.
 - Wireless and lightweight.
 - *AC power adapter.
 - ON/OFF switch.
 - Interchangeable heads (x8).
 - Weight: ~1 kg

POWER SUPPLY

110-240V 50/60Hz, 5A-26VDC-1A

Turning the device on:

Set the power switch (6) located on the bottom of the device to the ON position, then touch the touch-sensitive power button (5) to start the massager.

Adjusting massage intensity:

The device has 30 intensity levels. The current level is shown on the speed indicator (1).

“+” button (3) – increases the intensity.

“-” button (4) – decreases the intensity.

Changing the attachments:

Turn off the device before changing the attachment. Firmly push the selected attachment into the mounting hole until it is securely fitted.

Turning the device off:

After use, turn off the device using the touch-sensitive power button (5), then set the power switch (6) to the OFF position.

Charging the device:

Before first use, fully charge the battery. The first charging cycle may take up to 6 hours. The standard charging time is usually over 2 hours.

It is not recommended to fully discharge the battery. To maintain its service life, recharge the device regularly.

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